

Liquid Refilling Machine Production Line Sidolim 550

Brief Introduction

This automatic production line is consisted of the buffer turntable, automatic bottle washing machine, automatic linear liquid filling & capping machine, automatic screwing machine, automatic tabling machine. Automatically finish bottle washing, filling, plugging, capping procedure, also it can be consisted the perfect production line with dryer machine.

Advantage

Outer material is 304 stainless steel, and the material of filling parts is 316l stainless steel, it complies with GMP regulation and can be suitable for several bottle specs.

Stable running, low noise, precise filling

Low bottle breakage ratio, high qualified product ratio, low energy consumption

The whole line can single control or interlock control

Specifications

Filling volume: 50-500ml

Productivity: 40-60 bottle/min

Filling precision: <2%

Capping rate: >98%

Power supply: 380VAC/50Hz

Production line length: >7m

Infusion Solution Feeding Washing Filling
Plugging Capping and Labeling Production Line



www.sidolim.com



sale@sidolim.com

All products can be customized according to the customers' requirement and real bottle sample and liquid, and all products should be subjected to any change without notification. Please don't hesitate to contact us if have any questions.



Silver Double Inc.
Website: www.sidolim.com
Email: sale@sidolim.com

techservice@sidolim.com